

Raf1 (Phospho-Ser43) Rabbit mAb

Catalog No: #14232

Package Size: #14232-1 50ul #14232-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Raf1 (Phospho-Ser43) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2182
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	Phospho-Raf1 (S43) Antibody detects endogenous levels of total Phospho-Raf1 (S43)
Immunogen Description	A synthesized peptide derived from human Phospho-Raf1 (S43)
Conjugates	Unconjugated
Other Names	C-Raf; CRAF; Protein kinase raf 1; RAF1; Raf proto-oncogene serine/threonine-protein kinase; kinase Raf1;
Accession No.	Uniprot:P04049
Calculated MW	73 kDa
SDS-PAGE MW	65-75 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000

Images



All lanes : Phospho-c Raf (S43) Rabbit mAb at 1/1k dilution

Lane 1 : MCF7 whole cell lysates

Lane 2 : Hela whole cell lysates

Lane 3 : HepG2 whole cell lysates

Lane 4 : 293T whole cell lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 73 kDa

Observed band size: 65-75 kDa

Exposure time: 10 seconds



All lanes : Phospho-c Raf (S43) Rabbit mAb at 1/1k dilution

Lane 1 : 3T3 whole cell lysates

Lane 2 : Mouse Heart lysates

Lane 3 : Rat Brain lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 73 kDa

Observed band size: 65-75 kDa

Exposure time: 10 seconds



All lanes : Phospho-c Raf (S43) Rabbit mAb at 1/1k dilution

Lane 1 : MCF7 whole cell lysates

Lane 2 : MCF7 treated with 200nM PMA for 30min whole cell lysates

Lane 3 : MCF7 treated with Lambda Protein Phosphatase for 30min whole cell lysates

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution

Predicted band size: 73 kDa

Observed band size: 65-75 kDa

Exposure time: 7 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.