

beta Catenin (Phospho-Ser33/Ser37) Rabbit mAb

Catalog No: #14233

Package Size: #14233-1 50ul #14233-2 100ul


Description

Product Name	beta Catenin (Phospho-Ser33/Ser37) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2294
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Rat
Specificity	Phospho-beta Catenin (S33/S37) Antibody detects endogenous levels of total Phospho-beta Catenin (S33/S37)
Immunogen Description	A synthesized peptide derived from human Phospho-beta Catenin (S33/S37)
Conjugates	Unconjugated
Other Names	CATNB; CTNB1; CTNNB; catenin beta; beta catenin;
Accession No.	Uniprot:P35222
Calculated MW	Predicted band size: 85 kDa
SDS-PAGE MW	Observed band size: 92 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

Images



All lanes : beta Catenin (Phospho-Ser33/Ser37) Rabbit mAb at 1/1k dilution
Lane 1 :293T cell lysate
Lane 2 : 293T treated with 100nM Calyculin A for 30min whole cell lysates
lysate/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 85 kDa
Observed band size: 92 kDa
Exposure time: 9 seconds

kDa
140 —
115 —
80 —
65 —
50 —
40 —
30 —

1

beta Catenin (Phospho-Ser33/Ser37)

All lanes : beta Catenin (Phospho-Ser33/Ser37) Rabbit mAb at 1/1k dilution
Lane 1 : Rat brain lysate lysate/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 85 kDa
Observed band size: 92 kDa
Exposure time: 4 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.