

STAT3 (Phospho-Tyr705) Rabbit mAb

Catalog No: #14235

Package Size: #14235-1 50ul #14235-2 100ul



Description

Product Name	STAT3 (Phospho-Tyr705) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human Mouse Rat
Specificity	Phospho-STAT3 (Y705) Antibody detects endogenous levels of total Phospho-STAT3 (Y705)
Immunogen Description	A synthesized peptide derived from human STAT3
Conjugates	Unconjugated
Other Names	APRF; Stat3; HIES; Acute-phase response factor;
Accession No.	Uniprot:P40763
Calculated MW	88kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:5000~1:10000
IHC:1:50~1:100
ICC/IF:1:50~1:000
IP:1:30

Published Papers

Pei Gong;Pei Gong;Lijiao Zhao;Lijiao Zhao;Yunlong Ma;Yunlong Ma;Qiuting Shu;Qiuting Shu;Hui Sun;Hui Sun;Jing Lu;Jing Lu;Fanhua Meng;Fanhua Meng;Fang Wan;Fang Wan et al., AHR Agonist ITE Boosted PD1 Antibodyβ s Effects by Inhibiting Myeloid-Derived Cells Suppressive Cells in an Orthotopic Mouse Glioma Model, , (2025)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.