

Lyn (Phospho-Tyr396) Rabbit mAb

Catalog No: #14243

Package Size: #14243-1 50ul #14243-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Lyn (Phospho-Tyr396) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2011
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse
Specificity	Phospho-Lyn (Y396) Antibody detects endogenous levels of total Phospho-Lyn (Y396)
Immunogen Description	A synthesized peptide derived from human Phospho-Lyn (Y396)
Conjugates	Unconjugated
Other Names	JTK8; LYN; Tyrosine-protein kinase Lyn;
Accession No.	Uniprot:P07948
Calculated MW	Predicted band size: 59 kDa
SDS-PAGE MW	Observed band size: 56 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

Images



All lanes : Lyn (Phospho-Tyr396) Rabbit mAb at 1/1k dilution
 Lane 1 : Jurkat cell lysate
 Lane 2 : Jurkat cell lysate treated with Pervanadate.
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 59 kDa Observed band size: 56 kDa
 Exposure time: 4 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.