

IRF3 (Phospho-Ser386) Rabbit mAb

Catalog No: #14244

Package Size: #14244-1 50ul #14244-2 100ul


Description

Product Name	IRF3 (Phospho-Ser386) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	Phospho-IRF3 (S386) Antibody detects endogenous levels of total Phospho-IRF3 (S386)
Immunogen Description	A synthesized peptide derived from human Phospho-IRF3 (S386)
Conjugates	Unconjugated
Other Names	IRF3; Interferon regulatory factor 3;
Accession No.	Uniprot:Q14653
Calculated MW	47kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:1000~1:2000ICC/IF:1:50~1:200

Images



Western blot analysis of Phospho-IRF3 (S386) expression in MCF7 cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.