

MLKL (Phospho-Ser358) Rabbit mAb

Catalog No: #14251

Package Size: #14251-1 50ul #14251-2 100ul

Description


Product Name	MLKL (Phospho-Ser358) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1963
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human
Specificity	Phospho-MLKL (S358) Antibody detects endogenous levels of total Phospho-MLKL (S358)
Immunogen Description	A synthesized peptide derived from human Phospho-MLKL (S358)
Conjugates	Unconjugated
Other Names	hMLKL; Mlkl;
Accession No.	Uniprot:Q8NB16
Calculated MW	Predicted band size: 54 kDa
SDS-PAGE MW	Observed band size: 54 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : MLKL (Phospho-Ser358) Rabbit mAb at 1/1k dilution
 Lane 1 : HT-29 cell lysate
 Lane 2 : HT-29 cell lysate treated with 20 ng/ml TNF- α , 100 nM SM-164 and 20 μ M z-VAD
 Lysates/proteins at 20 μ g per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 54 kDa
 Observed band size: 54 kDa
 Exposure time: 9 seconds



Formalin-fixed, paraffin-embedded human small intestinal carcinoma tissue stained for MLKL (Phospho-Ser358) using 14251 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.