

KAP1 (Phospho-Ser824) Rabbit mAb

Catalog No: #14256

Package Size: #14256-1 50ul #14256-2 100ul

Description


Product Name	KAP1 (Phospho-Ser824) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1057
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human
Specificity	Phospho-KAP1 (S824) Antibody detects endogenous levels of total Phospho-KAP1 (S824)
Immunogen Description	A synthesized peptide derived from human KAP1
Conjugates	Unconjugated
Other Names	KAP 1; KRAB associated protein 1; KRIP1; RNF96; TF1B; TIF1 beta; TIF1B; Trim28; Tripartite motif containing 28;
Accession No.	Uniprot:Q13263
Calculated MW	Predicted band size: 89 kDa
SDS-PAGE MW	Observed band size: 110 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : KAP1 (Phospho-Ser824) Rabbit mAb at 1/1k dilution
Lane 1 :HeLa cell lysate
Lane 2 : HeLa treated with Etoposide cell lysate
Lysates/proteins at 20 µg per lane.
Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 89 kDa
Observed band size: 110 kDa
Exposure time: 6 seconds



Formalin-fixed, paraffin-embedded human tonsil tissue stained for KAP1 (Phospho-Ser824) using 14256 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.