

# Met (c-Met) (Phospho-Tyr1349) Rabbit mAb

Catalog No: #14257

Package Size: #14257-1 50ul #14257-2 100ul


## Description

|                       |                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------|
| Product Name          | Met (c-Met) (Phospho-Tyr1349) Rabbit mAb                                                            |
| Host Species          | Rabbit                                                                                              |
| Clonality             | Monoclonal                                                                                          |
| Isotype               | Rabbit IgG                                                                                          |
| Purification          | Affinity-chromatography                                                                             |
| Applications          | WB IHC IP                                                                                           |
| Species Reactivity    | Human                                                                                               |
| Specificity           | Phospho-Met (c-Met) (Y1349) Antibody detects endogenous levels of total Phospho-Met (c-Met) (Y1349) |
| Immunogen Description | A synthesized peptide derived from human c-Met                                                      |
| Conjugates            | Unconjugated                                                                                        |
| Other Names           | AUTS9; D249; HGF; HGF receptor; HGFR; MET; Par4; RCCP2;                                             |
| Accession No.         | Uniprot:P08581                                                                                      |
| Calculated MW         | 156kDa                                                                                              |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.  |
| Storage               | Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.  |

## Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

## Images



Western blot analysis of Phospho-Met (c-Met) (Y1349) expression in HeLa cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.