

IκB alpha (Phospho-Ser36) Rabbit mAb

Catalog No: #14260

Package Size: #14260-1 50ul #14260-2 100ul

Description


Product Name	IκB alpha (Phospho-Ser36) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1657
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human Mouse
Specificity	Phospho-IκB alpha (S36) Antibody detects endogenous levels of total Phospho-IκB alpha (S36)
Immunogen Description	A synthesized peptide derived from human Phospho-IκB alpha (S36)
Conjugates	Unconjugated
Other Names	I kappa B alpha; IkappaBalphα; IκB-α; IKBA; IκBalphα; MAD3; NF kappa B inhibitor alpha; NFKBI; NFKBIA;
Accession No.	Uniprot:P25963
Calculated MW	Predicted band size: 36 kDa
SDS-PAGE MW	Observed band size: 39 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : IκB alpha (Phospho-Ser36) Rabbit mAb at 1/1k dilution
Lane 1 : HeLa cell lysate treated with 100nM Calyculin A for 30min whole cell lysates
Lane 2 : HeLa cell lysate
Lysates/proteins at 20 μg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 36 kDa Observed band size: 39 kDa
Exposure time: 6 seconds



Formalin-fixed, paraffin-embedded human colon tissue stained for IKB alpha (Phospho-Ser36) using 14260 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.