

MLKL (Phospho-Ser345) Rabbit mAb

Catalog No: #14264

Package Size: #14264-1 50ul #14264-2 100ul

Description


Product Name	MLKL (Phospho-Ser345) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1962
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Mouse
Specificity	Phospho-MLKL (S345) Antibody detects endogenous levels of total Phospho-MLKL (S345)
Immunogen Description	A synthesized peptide derived from human Phospho-MLKL (S345)
Conjugates	Unconjugated
Other Names	hMLKL ; Mixed lineage kinase domain like; Mixed lineage kinase domain like protein; Mixed lineage kinase domain like pseudokinase;
Accession No.	Uniprot:Q9D2Y4
Calculated MW	Predicted band size: 54 kDa
SDS-PAGE MW	Observed band size: 54 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : MLKL (Phospho-Ser345) Rabbit mAb at 1/1k dilution
Lane 1 : L929 cell lysate
Lane 2 : L929 cell lysate treated with 20 ng/ml TNF-α, 100 nM SM-164 and 20 μM z-VAD
Lysates/proteins at 20 μg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 54 kDa
Observed band size: 54 kDa
Exposure time: 7 seconds



Formalin-fixed, paraffin-embedded mouse colon tissue stained for MLKL (Phospho-Ser345) using 14264 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.