

TAK1(Phospho-Ser439) Rabbit mAb

Catalog No: #14266

Package Size: #14266-1 50ul #14266-2 100ul


Description

Product Name	TAK1(Phospho-Ser439) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR2176
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-TAK1 (S439) Antibody detects endogenous levels of total Phospho-TAK1 (S439)
Immunogen Description	A synthesized peptide derived from human Phospho-TAK1 (S439)
Conjugates	Unconjugated
Other Names	MAP3K 7; MEKK7; Mitogen activated protein kinase kinase kinase 7; TAK1; TGF beta activated kinase 1; TGF1a;
Accession No.	Uniprot:O43318
Calculated MW	Predicted band size: 67 kDa
SDS-PAGE MW	Observed band size: 72 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000

Images



All lanes :TAK1(Phospho-Ser439) Rabbit mAb at 1/1k dilution
Lane 1 :K562 cell lysate
Lane 2 : THP-1 cell lysate
Lane 3 : Raji cell lysate
Lane 4 : MOLT-4 cell lysate
Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 67 kDa
Observed band size: 72 kDa
Exposure time: 13 seconds



All lanes :TAK1(Phospho-Ser439) Rabbit mAb at 1/1k dilution
 Lane 1 :HeLa treated with Calyculin A + IL-1 beta cell lysate
 Lane 2 : Untreated lysate
 Lysates/proteins at 20 µg per lane.
 SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 67 kDa Observed band size: 72 kDa
 Exposure time: 9 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.