

ACC(Phospho-Ser79) Rabbit mAb

Catalog No: #14267

Package Size: #14267-1 50ul #14267-2 100ul


Description

Product Name	ACC(Phospho-Ser79) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1413
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-ACC(S79) Antibody detects endogenous levels of total Phospho-ACC(S79)
Immunogen Description	A synthesized peptide derived from human Phospho-ACC(S79)
Conjugates	Unconjugated
Other Names	ACAC; ACACA; ACACB; ACC alpha; ACC; ACC beta; ACC1; ACC2; ACCA; ACCB; Acetyl CoA carboxylase 1; Acetyl CoA carboxylase 2; COA1; COA2;
Accession No.	Uniprot:Q13085
Calculated MW	Predicted band size: 266 kDa
SDS-PAGE MW	Observed band size: 265 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.


Application Details

WB: 1:500-1:2000

Images



All lanes :ACC(Phospho-Ser79) Rabbit mAb at 1/1k dilution
Lane 1 :A431 cell lysate
Lane 2 : A431 treated with Lambda phosphatase for 30min cell lysate
Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 266 kDa
Observed band size: 265 kDa
Exposure time: 8 seconds



All lanes :ACC(Phospho-Ser79) Rabbit mAb at 1/1k dilution
Lane 1 :293 cell lysate
Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 266 kDa
Observed band size: 265 kDa
Exposure time: 5 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.