

CBL(Phospho-Ser669) Rabbit mAb

Catalog No: #14272

Package Size: #14272-1 50ul #14272-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CBL(Phospho-Ser669) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1097
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-CBL (S669) Antibody detects endogenous levels of total Phospho-CBL (S669)
Immunogen Description	A synthesized peptide derived from human Phospho-CBL (S669)
Conjugates	Unconjugated
Other Names	Casitas B lineage lymphoma proto oncogene; cbl; CBL2; E3 ubiquitin protein ligase CBL; Oncogene CBL2; Proto oncogene c CBL; RING finger protein 55; RNF55; Signal transduction protein CBL;
Accession No.	Uniprot:P22681
Calculated MW	Predicted band size: 100 kDa
SDS-PAGE MW	Observed band size: 120 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB: 1:500-1:2000

Images



All lanes :CBL(Phospho-Ser669) Rabbit mAb at 1/1k dilution
 Lane 1 :HeLa cell lysate
 Lane 2 : HeLa treated with pervanadate lysate
 Lysates/proteins at 20 μ g per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 100 kDa
 Observed band size: 120 kDa
 Exposure time: 9 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.