

Nucleolin (Phospho-Thr84) Rabbit mAb

Catalog No: #14279

Package Size: #14279-1 50ul #14279-2 100ul



SAB[®]
Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Nucleolin (Phospho-Thr84) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1020
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human
Specificity	Phospho-Nucleolin (T84) Antibody detects endogenous levels of total Phospho-Nucleolin (T84)
Immunogen Description	A synthesized peptide derived from human Phospho-Nucleolin (T84)
Conjugates	Unconjugated
Other Names	NCL; Nucl; Nucleolin; Protein C23;
Accession No.	Uniprot:P19338
Calculated MW	Predicted band size: 77 kDa
SDS-PAGE MW	Observed band size: 100 kDa
Concentration	2mg/ml
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes :Nucleolin (Phospho-Thr84) Rabbit mAb at 1/1k dilution
Lane 1 :K562 cell lysate
Lane 2 : K562 treated with 100nM Calyculin A for 30min whole cell lysates
Lysates/proteins at 20 µg per lane.
SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 77 kDa Observed band size: 100 kDa
Exposure time: 5 seconds



Formalin-fixed, paraffin-embedded human tonsil tissue stained for Nucleolin (Phospho-Thr84) using 14279 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.