

JAK1 (Phospho-Tyr1034/Tyr1035) Rabbit mAb

Catalog No: #14286

Package Size: #14286-1 50ul #14286-2 100ul

Description


Product Name	JAK1 (Phospho-Tyr1034/Tyr1035) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR1138
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,IHC
Species Reactivity	Human
Specificity	Phospho-JAK1 (Y1034 + Y1035) Antibody detects endogenous levels of total Phospho-JAK1 (Y1034 + Y1035)
Immunogen Description	A synthesized peptide derived from human Phospho-JAK1 (Y1034 + Y1035)
Conjugates	Unconjugated
Other Names	JAK 1; JAK 1A; JAK 1B; JAK1; JAK1A; JAK1B; JTK3;
Accession No.	Uniprot:P23458
Calculated MW	Predicted band size: 133 kDa
SDS-PAGE MW	Observed band size: 133 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200

Images



All lanes : JAK1 (Phospho-Tyr1034/Tyr1035) Rabbit mAb at 1/1k dilution Lane 1 :Ramos treated with pervanadate cell lysate.Lysates/proteins at 20 µg per lane.SecondaryAll lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilutionPredicted band size: 133 kDa Observed band size: 133 kDaExposure time: 5 seconds



Formalin-fixed, paraffin-embedded human gastric adenocarcinoma tissue stained for JAK1 (Phospho-Tyr1034/Tyr1035) using 14286 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.