

Tau (Phospho-Ser202) Rabbit mAb

Catalog No: #14290

Package Size: #14290-1 50ul #14290-2 100ul


Description

Product Name	Tau (Phospho-Ser202) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human Mouse Rat
Specificity	Phospho-Tau (S202) Antibody detects endogenous levels of total Phospho-Tau (S202)
Immunogen Description	A synthesized peptide derived from human Phospho-Tau (S202)
Conjugates	Unconjugated
Other Names	MAPT; Microtubule-associated protein tau; MTBT1; Neurofibrillary tangle protein; Paired helical filament-tau; PHF-tau
Accession No.	Uniprot:P10636
Calculated MW	50-70kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of Phospho-Tau (S202) expression in mouse hippocampus cell lysate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.