

Histone H3 (Phospho Thr3) Rabbit mAb

Catalog No: #14313

Package Size: #14313-1 50ul #14313-2 100ul

Description

Product Name	Histone H3 (Phospho Thr3) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Purification	Protein A
Applications	WB;IHC;IF;IP;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Phospho-Histone H3 (T3) Antibody detects endogenous levels of Histone H3 protein only when phosphorylated at T3.
Target Name	Histone H3
Modification	Phospho
Other Names	HIST1H3A;H3FA;HIST1H3B;H3FL;HIST1H3C;H3FC;HIST1H3D;H3FB;HIST1H3E;H3FD;HIST1H3F;H3FI;HIST1H3G;H3FH;HIST1H3H;H3FK;HIST1H3I;H3FF;HIST1H3J;H3FJ;Histone H3.1;Histone H3/a;Histone H3/b;Histone H3/c;Histone H3/d;Histone H3
Accession No.	P68431
Calculated MW	15kDa
SDS-PAGE MW	15kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details

IHC 1:200-1:1000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.