

MEK1/2 (Phospho Ser217/221) Rabbit mAb

Catalog No: #14315

Package Size: #14315-1 50ul #14315-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEK1/2 (Phospho Ser217/221) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Purification	Protein A
Applications	WB;IHC;IF;IP;ELISA
Species Reactivity	Human;Mouse;Rat
Target Name	MEK1/2
Modification	Phospho
Accession No.	P36507;Q02750
Calculated MW	44kDa
SDS-PAGE MW	44kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details

IHC 1:1000-1:5000;WB 1:2000-1:10000;IF 1:200-1:1000;ELISA 1:5000-1:20000;IP 1:50-1:200;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.