

Syk (Phospho-Tyr525/526) Rabbit mAb

Catalog No: #14325

Package Size: #14325-1 50ul #14325-2 100ul


Description

Product Name	Syk (Phospho-Tyr525/526) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR4365
Purification	Protein A
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	Antibody detects endogenous levels of syk protein only when phosphorylated at Tyr525/526 of human Syk or Tyr519/520 of mouse Syk.
Target Name	Syk
Modification	Phospho
Other Names	SYK;Tyrosine-protein kinase SYK;Spleen tyrosine kinase;p72-Syk
Accession No.	P43405
Calculated MW	Predicted band size: 72 kDa
SDS-PAGE MW	Observed band size: 72 kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details

WB: 1:500-1:2000

Images



All lanes : Syk (Phospho-Tyr525/526) Rabbit mAb at 1/1k dilution
Lane 1 : Ramos whole cell lysates
Lane 2 : Ramos treated with 12 ug/ml anti-IgM for 2 minutes whole cell lysates
Lane 3 : C6 whole cell lysates
Lane 4 : C6 treated with 10 nM Na3VO4 for 1 h whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 72 kDa
Observed band size: 72 kDa
Exposure time: 5 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.