

Nucleophosmin Rabbit mAb

Catalog No: #14329

Package Size: #14329-1 50ul #14329-2 100ul



Signalway Antibody

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	Nucleophosmin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Purification	Protein A
Applications	WB,IHC
Species Reactivity	Human;Mouse;Rat
Specificity	Endogenous
Target Name	Nucleophosmin
Modification	Unmodified
Other Names	NPM1;NPM;Nucleophosmin;NPM;Nucleolar phosphoprotein B23;Nucleolar protein NO38;Numatrin
Accession No.	P06748
Calculated MW	Predicted band size: 33 kDa
SDS-PAGE MW	Observed band size: 38 kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details


WB: 1:500-1:2000

IHC: 1:50-1:200


Images



All lanes : Nucleophosmin Rabbit mAb at 1/1k dilution
Lane 1 : Hela whole cell lysates
Lane 2 : HepG2 whole cell lysates
Lane 3 : HCT-116 whole cell lysates
Lane 4 : C2C12 whole cell lysates
Lane 5 : Raw264.7 whole cell lysates
Lysates/proteins at 20 µg per lane.
Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
Predicted band size: 33 kDa
Observed band size: 38 kDa
Exposure time: 4 seconds



Formalin-fixed, paraffin-embedded human liver tissue stained for Nucleophosmin using 14329 at 1/100 dilution in immunohistochemical analysis.



Formalin-fixed, paraffin-embedded mouse liver tissue stained for Nucleophosmin using 14329 at 1/100 dilution in immunohistochemical analysis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.