

PERK (Phospho-Thr980) Rabbit mAb

Catalog No: #14330

Package Size: #14330-1 50ul #14330-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PERK (Phospho-Thr980) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR4369
Purification	Protein A
Applications	WB
Species Reactivity	Mouse;Rat
Specificity	PERK (Phospho Thr980) Antibody detects endogenous levels of PERK protein only when phosphorylated at T980.
Target Name	PERK
Modification	Phospho
Other Names	EIF2AK3;PEK;PERK;Eukaryotic translation initiation factor 2-alpha kinase 3;PRKR-like endoplasmic reticulum kinase;Pancreatic eIF2-alpha kinase;HsPEK
Accession No.	Q9NZJ5
Calculated MW	Predicted band size: 125 kDa
SDS-PAGE MW	Observed band size: 125 kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details

WB: 1:500-1:2000

Images



All lanes : PERK (Phospho-Thr980) Rabbit mAb at 1/1k dilution
 Lane 1 : Rat pancreas lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 125 kDa
 Observed band size: 125 kDa
 Exposure time: 6 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.