

IRF3 (Phospho-Ser396) Rabbit mAb

Catalog No: #14331

Package Size: #14331-1 50ul #14331-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IRF3 (Phospho-Ser396) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR4370
Purification	Protein A
Applications	WB
Species Reactivity	Human;Mouse
Specificity	IRF3 (Phospho Ser396) Antibody detects endogenous levels of IRF3 protein only when phosphorylated at S396.
Target Name	IRF-3
Modification	Phospho
Other Names	IRF3;Interferon regulatory factor 3;IRF-3
Accession No.	Q14653
Calculated MW	Predicted band size: 47 kDa
SDS-PAGE MW	Observed band size: 60 kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details

WB: 1:500-1:2000

Images



All lanes : IRF3 (Phospho-Ser396) Rabbit mAb at 1/1k dilution
 Lane 1 : HT29 whole cell lysates
 Lane 2 : HT29 treated with Lambda Protein Phosphatase for 30min whole cell
 Lane 3 : 3T3 whole cell lysates
 Lysates/proteins at 20 µg per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 47 kDa Observed band size: 60 kDa
 Exposure time: 9 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.