

PLC γ 1 (Phospho-Tyr783) Rabbit mAb

Catalog No: #14335

Package Size: #14335-1 50ul #14335-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PLC γ 1 (Phospho-Tyr783) Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Clone No.	SR4374
Purification	Protein A
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	PLC γ 1 (Phospho Tyr783) Monoclonal Antibody detects endogenous levels of PLC γ 1 protein only when phosphorylated at Y783.
Target Name	PLCG1
Modification	Phospho
Other Names	PLCG1;PLC1;1-phosphatidylinositol 4;5-bisphosphate phosphodiesterase gamma-1;PLC-148;Phosphoinositide phospholipase C-gamma-1;Phospholipase C-II;PLC-II;Phospholipase C-gamma-1;PLC-gamma-1
Accession No.	P19174
Calculated MW	Predicted band size: 149 kDa
SDS-PAGE MW	Observed band size: 149 kDa
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Storage	Store at -20°C /1 year

Application Details

WB: 1:500-1:2000

Images



All lanes : PLC γ 1 (Phospho-Tyr783) Rabbit mAb at 1/1k dilution
 Lane 1 : HeLa whole cell lysates
 Lane 2 : HeLa treated with 100ng/ml PDGF for 5 minutes whole cell lysates
 Lane 3 : 3T3 whole cell lysates
 Lane 4 : 3T3 treated with 100ng/ml PDGF for 5 minutes whole cell lysates
 Lane 5 : C6 whole cell lysates
 Lane 6 : C6 treated with 100ng/ml PDGF for 5 minutes whole cell lysates
 Lysates/proteins at 20 μ g per lane.
 Secondary All lanes : Goat Anti-Rabbit IgG H&L (HRP) at 1/20000 dilution
 Predicted band size: 149 kDa Observed band size: 149 kDa
 Exposure time: 11 seconds

Note: This product is for in vitro research use only and is not intended for use in humans or animals.