

HDAC4(Ab-632) Antibody

Catalog No: #21141

Package Size: #21141-1 50ul #21141-2 100ul

Description

Product Name	HDAC4(Ab-632) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total HDAC4 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa. 630~634 (A-Q-S-S-P) derived from Human HDAC4.
Conjugates	Unconjugated
Target Name	HDAC4
Other Names	HD4; Histone deacetylase 4;
Accession No.	Swiss-Prot: P56524NCBI Protein: NP_006028.2
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 140kd


Western blotting: 1:500~1:1000

Immunofluorescence: 1:100~1:200


Images



Western blot analysis of extracts from HepG2 cells using HDAC4(Ab-632) Antibody #21141.



Immunofluorescence staining of methanol-fixed HeLa cells using HDAC4(Ab-632) Antibody #21141.



Western blotting analysis using HDAC4(Ab-632) Antibody #21141.

Background

Responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes. Involved in muscle maturation via its interaction with the myocyte enhancer factors such as MEF2A, MEF2C and MEF2D.

Wang AH, et al. (2000) Mol Cell Biol. 20(18): 6904-6912.

Grozinger CM, et al. (2000) Proc Natl Acad Sci U S A. 97(14): 7835-7840.

Published Papers

el at., Anticancer effects of sodium butyrate on hepatocellular carcinoma cells in vitro. In Int J Mol Med on 2013 Apr by Hong-Gang Wang, Xiao-Dan Huang, et al..PMID:

23440283, , (2013)

[PMID:23440283](#)

el at., Expression and localization of calmodulin-related proteins in brain, heart and kidney from spontaneously hypertensive rats.In Biochem Biophys Res Commun on 2016 Jan 15 by Satoshi Kameshima, Muneyoshi Okada et al..PMID:26697749, , (2016)

[PMID:26697749](#)

el at., Thrombospondin-1 might be a therapeutic target to suppress RB cells by regulating the DNA double-strand breaks repair. In Oncotarget on 2016 Feb 2 by Pei Chen, Na Yu, et al..PMID:26756218, , (2016)

[PMID:26756218](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.