

b-Enolase(ENO-3) Antibody

Catalog No: #21656

Package Size: #21656-1 50ul #21656-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	b-Enolase(ENO-3) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous level of total b-Enolase(ENO-3) protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.49~53(L-R-D-G-D) derived from Human b-Enolase(ENO-3).
Conjugates	Unconjugated
Target Name	b-Enolase(ENO-3)
Other Names	MSE; GSD13;
Accession No.	Swiss-Prot: P13929NCBI Protein: NP_001967.3
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 45kd

Western blotting: 1:500~1:1000

Images



Western blot analysis of extracts from 3T3, HepG2 and C2C12 cells using b-Enolase(ENO-3) Antibody #21656.



Western blot analysis of extracts from
HepG2oO JKoO HeLa,3T3, C2C12 and 293T cells using
b-Enolase(ENO-3) Antibody

Background

Appears to have a function in striated muscle development and regeneration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.