

S6 Ribosomal Protein (Ab-235/236) antibody

Catalog No: #21682

Package Size: #21682-1 50ul #21682-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	S6 Ribosomal Protein (Ab-235/236) antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total S6 Ribosomal protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa.233-238(R-L-S-S-L-R) derived from Human S6 Ribosomal Protein.
Conjugates	Unconjugated
Target Name	S6 Ribosomal Protein
Other Names	NP33, RPS6, RS6,
Accession No.	Swiss-Prot:P62753Gene ID:6194
Concentration	1.0mg/ml

Application Details

Predicted MW: 32kd

Western blotting: 1: 500~1: 1000

Images



Western blot analysis of extracts from HeLa, 293 and PC12 cells using S6 Ribosomal Protein (Ab-235/236) Antibody #21682.

Background

May play an important role in controlling cell growth and proliferation through the selective translation of particular classes of mRNA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.