

## CYFIP2 antibody

Catalog No: #22012

Package Size: #22012 100ul

## Description

|                       |                                                                                                                    |
|-----------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name          | CYFIP2 antibody                                                                                                    |
| Host Species          | Rabbit                                                                                                             |
| Clonality             | Polyclonal                                                                                                         |
| Purification          | Purified by antigen-affinity chromatography.                                                                       |
| Applications          | WB IHC                                                                                                             |
| Species Reactivity    | Hu                                                                                                                 |
| Immunogen Type        | Recombinant protein                                                                                                |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 346 and 610 of CYFIP2 |
| Conjugates            | Unconjugated                                                                                                       |
| Target Name           | CYFIP2                                                                                                             |
| Accession No.         | Swiss-Prot:Q96F07Gene ID:26999                                                                                     |
| Concentration         | 1mg/ml                                                                                                             |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.     |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                          |


## Application Details

Predicted MW: 146kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

## Images



Sample (30 ug of whole cell lysate)  
A: Molt-4  
5% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded NCIN87 xenograft, using CYFIP2 antibody at 1: 100 dilution.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.