

DLST antibody

Catalog No: #22024

Package Size: #22024 100ul


Description

|                       |                                                                                                                        |
|-----------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name          | DLST antibody                                                                                                          |
| Host Species          | Rabbit                                                                                                                 |
| Clonality             | Polyclonal                                                                                                             |
| Purification          | Purified by antigen-affinity chromatography.                                                                           |
| Applications          | WB IHC IF                                                                                                              |
| Species Reactivity    | Hu                                                                                                                     |
| Immunogen Type        | Recombinant protein                                                                                                    |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 188 and 423 of human DLST |
| Conjugates            | Unconjugated                                                                                                           |
| Target Name           | DLST                                                                                                                   |
| Accession No.         | Swiss-Prot:P36957Gene ID:1743                                                                                          |
| Concentration         | 0.7mg/ml                                                                                                               |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.         |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                              |


Application Details

|                                  |
|----------------------------------|
| Predicted MW: 49kd               |
| Western blotting: 1:500-1:3000   |
| Immunohistochemistry: 1:50-1:500 |
| Immunofluorescence: 1:100-1:200  |


Images



Sample (30 ug of whole cell lysate)  
A: A431  
10% SDS PAGE  
Primary antibody diluted at 1: 1000



A:Immunohistochemical analysis of paraffin-embedded RT2 xenograft, using DLST antibody at 1: 500 dilution.  
B:Immunohistochemical analysis of paraffin-embedded Cal27 xenograft, using DLST antibody at 1: 500 dilution.



Costained with Hoechst 33342

Immunofluorescence analysis of methanol-fixed A431, using DLST antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.