

IFIT3 antibody

Catalog No: #22040

Package Size: #22040 100ul

Description

Product Name	IFIT3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 230 and 470 of Human IFIT3
Conjugates	Unconjugated
Target Name	IFIT3
Accession No.	Swiss-Prot:O14879Gene ID:3437
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 56kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images



Sample (30 ug of whole cell lysate)
A: 293T
7.5% SDS PAGE
Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded Hela xenograft, using IFIT3 antibody at 1: 500 dilution.

Published Papers

el at., NF- κ B-Dependent IFIT3 Induction by HBx Promotes Hepatitis B Virus Replication. In Front Microbiol on 2019 Oct 11 by Xu F, Song H, et al..PMID:31681236, , (2019)

[PMID:31681236](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.