

GnT-III antibody

Catalog No: #22062

Package Size: #22062 100ul

Description

Product Name	GnT-III antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 260 and 495 of Human MGAT3
Conjugates	Unconjugated
Target Name	GnT-III
Accession No.	Swiss-Prot:Q09327Gene ID:4248
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 61kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:50-1:500

Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)
A: Hep G2
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SW480 xenograft, using MGAT3 antibody at 1: 500 dilution.



Immunofluorescence analysis of methanol-fixed HeLa, using GnT-III antibody at 1: 200 dilution.

Background

There are believed to be over 100 different glycosyltransferases involved in the synthesis of protein-bound and lipid-bound oligosaccharides. The enzyme encoded by this gene transfers a GlcNAc residue to the beta-linked mannose of the trimannosyl core of N-linked oligosaccharides and produces a bisecting GlcNAc. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.