

## KRR1 antibody

Catalog No: #22069

Package Size: #22069 100ul

## Description

|                       |                                                                                                                 |
|-----------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name          | KRR1 antibody                                                                                                   |
| Host Species          | Rabbit                                                                                                          |
| Clonality             | Polyclonal                                                                                                      |
| Purification          | Purified by antigen-affinity chromatography.                                                                    |
| Applications          | WB IHC IF                                                                                                       |
| Species Reactivity    | Hu                                                                                                              |
| Immunogen Type        | Recombinant protein                                                                                             |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 47 and 314 of KRR1 |
| Conjugates            | Unconjugated                                                                                                    |
| Target Name           | KRR1                                                                                                            |
| Accession No.         | Swiss-Prot:Q13601Gene ID:11103                                                                                  |
| Concentration         | 1mg/ml                                                                                                          |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.  |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                       |

## Application Details


Predicted MW: 44kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:100-1:500

Immunofluorescence: 1:100-1:200


## Images




Sample (30 ug of whole cell lysate)  
A: H1299  
B: HeLa  
10% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SAS xenograft, using KRR1 antibody at 1: 500 dilution.



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using KRR1 antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.