

RAB6A antibody

Catalog No: #22075

Package Size: #22075 100ul

Description

Product Name	RAB6A antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 178 of RAB6A
Conjugates	Unconjugated
Target Name	RAB6A
Accession No.	Swiss-Prot:P20340Gene ID:5870
Concentration	0.21 mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 24kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:100-1:500

Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)
A: H1299
12% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded MDA-MB-468 xenograft, using RAB6A antibody at 1: 500 dilution.



Confocal immunofluorescence analysis (Olympus FV10i) of methanol-fixed HeLa, using RAB6A antibody (Green) at 1: 1000 dilution and alpha-tubulin antibody (Red) at 1: 2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.