

ZNF707 antibody

Catalog No: #22104

Package Size: #22104 100ul

Description

Product Name	ZNF707 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 121 and 345 of ZNF707
Conjugates	Unconjugated
Target Name	ZNF707
Accession No.	Swiss-Prot:Q96C28Gene ID:286075
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 43kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

Images



Sample (30 ug of whole cell lysate)
A: Hep G2
10% SDS PAGE
Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded H3347 xenograft, using ZNF707 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.