

Nectin 2 antibody

Catalog No: #22110

Package Size: #22110 100ul

Description


Product Name	Nectin 2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 270 and 509 of Human PVRL2
Conjugates	Unconjugated
Target Name	Nectin 2
Accession No.	Swiss-Prot:Q92692Gene ID:5819
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 58kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: 293T
10% SDS PAGE
Primary antibody diluted at 1: 500

Background

This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and

pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.