

CKLFSF5 antibody

Catalog No: #22137

Package Size: #22137 100ul

Description


Product Name	CKLFSF5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified by Protein A.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 2 and 15 of CKLFSF5
Conjugates	Unconjugated
Target Name	CKLFSF5
Other Names	CKLFSF5; FLJ37521
Accession No.	Swiss-Prot:Q96DZ9Gene ID:116173
Concentration	5.8mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 17kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: Raji
12% SDS PAGE
CKLFSF5 antibody diluted at 1: 500

Background

This gene belongs to the chemokine-like factor gene superfamily. This family of genes encodes multi-pass membrane proteins that are similar to both the chemokine and the transmembrane 4 superfamilies of signaling molecules. Alternate transcriptional splice variants of this gene, encoding different

isoforms, have been characterized. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.