

PECR antibody

Catalog No: #22145

Package Size: #22145 100ul

Description

Product Name	PECR antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 241 and 302 of human PECR
Conjugates	Unconjugated
Target Name	PECR
Accession No.	Swiss-Prot:Q9BY49Gene ID:55825
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 33kd


Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)
A: HeLa
B: Hep G2
12% SDS PAGE
Primary antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed A431, using PECC antibody at 1: 200 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.