

MANBA antibody

Catalog No: #22166

Package Size: #22166 100ul

Description


Product Name	MANBA antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 380 and 624 of Human MANBA
Conjugates	Unconjugated
Target Name	MANBA
Accession No.	Swiss-Prot:O00462Gene ID:4126
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 101kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: Hep G2
7.5% SDS PAGE
Primary antibody diluted at 1: 500

Background

This gene encodes a member of the glycosyl hydrolase 2 family. The encoded protein localizes to the lysosome where it is the final exoglycosidase in the pathway for N-linked glycoprotein oligosaccharide catabolism. Mutations in this gene are associated with beta-mannosidosis, a lysosomal storage

disease that has a wide spectrum of neurological involvement. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.