

SUGLG2 antibody

Catalog No: #22171

Package Size: #22171 100ul

Description

Product Name	SUGLG2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 204 of Human SUGLG2
Conjugates	Unconjugated
Target Name	SUGLG2
Accession No.	Swiss-Prot:Q96I99Gene ID:8801
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 47kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images



Sample(30ug whole cell lysate)

A: 293T

B: A431

C: H1299

D: HeLa S3


E: Hep G2

F: MOLT4


G: Raji

10% SDS PAGE

Primary antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed HeLa, using SUCLG2 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.