

MSTO1 antibody

Catalog No: #22173

Package Size: #22173 100ul

Description

Product Name	MSTO1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 34 and 297 of MSTO1
Conjugates	Unconjugated
Target Name	MSTO1
Other Names	DKFZp686B1757; DKFZp686I01261; FLJ10504; LST005; MST; RP11-29H23.3
Accession No.	Swiss-Prot:Q9BUK6Gene ID:55154
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 62kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample (30 ug of whole cell lysate)
A: H1299
7.5% SDS PAGE
MSTO1 antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded OVCA xenograft, using MSTO1 antibody at 1: 100 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.