

## FLJ23560 antibody

Catalog No: #22176

Package Size: #22176 100ul

## Description

|                       |                                                                                                                         |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name          | FLJ23560 antibody                                                                                                       |
| Host Species          | Rabbit                                                                                                                  |
| Clonality             | Polyclonal                                                                                                              |
| Purification          | Purified by antigen-affinity chromatography.                                                                            |
| Applications          | WB IHC                                                                                                                  |
| Species Reactivity    | Hu                                                                                                                      |
| Immunogen Type        | Recombinant protein                                                                                                     |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 182 and 413 of Human BBS10 |
| Conjugates            | Unconjugated                                                                                                            |
| Target Name           | FLJ23560                                                                                                                |
| Accession No.         | Swiss-Prot:Q8TAM1Gene ID:79738                                                                                          |
| Concentration         | 1mg/ml                                                                                                                  |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.          |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                               |


## Application Details

Predicted MW: 81kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

## Images



Sample (30 ug of whole cell lysate)  
A: Hep G2  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded NCI-N87 xenograft, using BBS10 antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.