

## DKFZP686A01247 antibody

Catalog No: #22177

Package Size: #22177 100ul

## Description

|                       |                                                                                                                           |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name          | DKFZP686A01247 antibody                                                                                                   |
| Host Species          | Rabbit                                                                                                                    |
| Clonality             | Polyclonal                                                                                                                |
| Purification          | Purified by antigen-affinity chromatography.                                                                              |
| Applications          | WB IHC                                                                                                                    |
| Species Reactivity    | Hu                                                                                                                        |
| Immunogen Type        | Recombinant protein                                                                                                       |
| Immunogen Description | Recombinant protein fragment contain a sequence corresponding to a region within amino acids 734 and 1057 of Human LIMCH1 |
| Conjugates            | Unconjugated                                                                                                              |
| Target Name           | DKFZP686A01247                                                                                                            |
| Accession No.         | Swiss-Prot:Q9UPQ0Gene ID:22998                                                                                            |
| Concentration         | 1mg/ml                                                                                                                    |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.            |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                                 |


## Application Details

Predicted MW: 119kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

## Images



Sample (30ug whole cell lysate)  
A: H1299  
7.5% SDS PAGE  
Primary antibody diluted at 1: 3000



Immunohistochemical analysis of paraffin-embedded A549 xenograft, using LIMCH1 antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.