

muscleblind-like 3 isoform G antibody

Catalog No: #22178

Package Size: #22178 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	muscleblind-like 3 isoform G antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 267 of Human MBNL3
Conjugates	Unconjugated
Target Name	muscleblind-like 3 isoform G
Accession No.	Swiss-Prot:Q9NUK0Gene ID:55796
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 39kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200


Images



Sample(30 ug whole cell lysate)
 A: A431
 B: Raji
 10% SDS PAGE
 Primary antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of paraformaldehyde-fixed A431, using MBNL3 antibody at 1: 200 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.