

OSGEP antibody

Catalog No: #22180

Package Size: #22180 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OSGEP antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 140 and 332 of Human OSGEP
Conjugates	Unconjugated
Target Name	OSGEP
Other Names	FLJ20411; GCPL1; KAE1; OSGEP1; PRSMG1
Accession No.	Swiss-Prot:Q9NPF4 Gene ID:55644
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 36kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images




Sample(30 ug whole cell lysate)

A: Raji

10% SDS PAGE

Primary antibody diluted at 1: 2000



Immunohistochemical analysis of paraffin-embedded Colon ca, using OSGEP antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.