

MMACHC antibody

Catalog No: #22182

Package Size: #22182 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MMACHC antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 188 of Human MMACHC
Conjugates	Unconjugated
Target Name	MMACHC
Accession No.	Swiss-Prot:Q9Y4U1Gene ID:25974
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 32kd

Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Immunofluorescence: 1:100-1:200

Images




Sample(30 ug whole cell lysate)

A: 293T


B: Raji

12% SDS PAGE

Primary antibody diluted at 1: 1000



A:Immunohistochemical analysis of paraffin-embedded C2C12 xenograft, using MMACHC antibody at 1: 500 dilution.
B:Immunohistochemical analysis of paraffin-embedded Hepatoma, using MMACHC antibody at 1: 500 dilution.



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using MMACHC antibody at 1: 100 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.