

OSBPL6 antibody

Catalog No: #22195

Package Size: #22195 100ul

Description

Product Name	OSBPL6 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide contain a sequence corresponding to a region within amino acids 872 and 934 of human OSBPL6
Conjugates	Unconjugated
Target Name	OSBPL6
Accession No.	Swiss-Prot:Q9BZF3 Gene ID:114880
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 106kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:100-1:250

Immunofluorescence: 1:100-1:200


Images




Sample (30 ug of whole cell lysate)
A: A431
7.5% SDS PAGE
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded SCC15 xenograft, using OSBPL6 antibody at 1: 100 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of paraformaldehyde-fixed HeLa, using OSBPL6 antibody at 1: 200 dilution.

Background

This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. Transcript variants encoding different isoforms have been identified. [provided by RefSeq]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.