

CGI-121 antibody

Catalog No: #22213

Package Size: #22213 100ul

Description

Product Name	CGI-121 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 159 of Human TPRKB
Conjugates	Unconjugated
Target Name	CGI-121
Accession No.	Swiss-Prot:Q9Y3C4Gene ID:51002
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 20kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:250

Images



Sample(30 ug whole cell lysate)
A: 293T
12% SDS PAGE
Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded Ca922 xenograft, using TPRKB antibody at 1: 250 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.