

## EVC2 antibody

Catalog No: #22218

Package Size: #22218 100ul

## Description

|                       |                                                                                                                |
|-----------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name          | EVC2 antibody                                                                                                  |
| Host Species          | Rabbit                                                                                                         |
| Clonality             | Polyclonal                                                                                                     |
| Purification          | Purified by antigen-affinity chromatography.                                                                   |
| Applications          | WB IF                                                                                                          |
| Species Reactivity    | Hu                                                                                                             |
| Immunogen Type        | Peptide                                                                                                        |
| Immunogen Description | Synthetic peptide contain a sequence corresponding to a region within amino acids 1246 and 1308 of EVC2        |
| Conjugates            | Unconjugated                                                                                                   |
| Target Name           | EVC2                                                                                                           |
| Accession No.         | Swiss-Prot:Q86UK5Gene ID:132884                                                                                |
| Concentration         | 1mg/ml                                                                                                         |
| Formulation           | Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative. |
| Storage               | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.                      |

## Application Details


Predicted MW: 148kd

Western blotting: 1:500-1:3000


Immunohistochemistry: 1:100-1:250

Immunofluorescence: 1:100-1:200


## Images



Sample (30 ug of whole cell lysate)  
A: H1299  
5% SDS PAGE  
EVC2 antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed HeLa, using EVC2 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.