

CCDC83 antibody

Catalog No: #22222

Package Size: #22222 100ul

Description

Product Name	CCDC83 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 163 and 407 of CCDC83
Conjugates	Unconjugated
Target Name	CCDC83
Accession No.	Swiss-Prot:Q8IWF9Gene ID:220047
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 52kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

Images



Sample (30 ug of whole cell lysate)
A: JurKat
10% SDS PAGE
CCDC83 antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Normal Colon, using CCDC83 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.