

TXNDC3 antibody

Catalog No: #22223

Package Size: #22223 100ul

Description


Product Name	TXNDC3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 270 and 520 of Human TXNDC3
Conjugates	Unconjugated
Target Name	TXNDC3
Accession No.	Swiss-Prot:Q8N427Gene ID:51314
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 67kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: Raji

7.5% SDS PAGE

Primary antibody diluted at 1: 1000

Background

TXNDC3 encodes a group II thioredoxin protein, composed of a thioredoxin domain and 3 NDP kinase domains, with testis-specific expression.[supplied by OMIM]

Note: This product is for in vitro research use only and is not intended for use in humans or animals.