

IL-17 receptor D antibody

Catalog No: #22228

Package Size: #22228 100ul

Description


Product Name	IL-17 receptor D antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 286 and 501 of Human IL17RD
Conjugates	Unconjugated
Target Name	IL-17 receptor D
Accession No.	Swiss-Prot:Q8NFM7Gene ID:54756
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 82kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)
A: 293T
7.5% SDS PAGE
Primary antibody diluted at 1: 500

Background

Fibroblast growth factors (FGFs; see MIM 603726) are secreted proteins involved in cellular proliferation, migration, differentiation, and survival. FGF activity is negatively regulated by members of the 'sprouty' family (e.g., SPRY1, MIM 602465). The SEF protein is a modulator of FGF

Note: This product is for in vitro research use only and is not intended for use in humans or animals.